

Product Datasheet

Recombinant Mouse FGF7 Protein orb527082

Reactivity	Mouse
Form/Appearance	Lyophilized
Storage	Lyophilized recombinant mouse FGF7 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rmFGF7 remains stable up to 2 weeks at 4°C or up to 3 months at -20°C.
Note	For research use only.
Protein Sequence	CNDMS PEQTA TSVNC SSPER HTRSY DYMEG GDIRV RRLFC RTQWY LRIDK RGKVK GTQEM KNSYN IMEIR TVAVG IVAIK GVESE YYLAM NKEGK LYAKK ECNED CNFKE LILEN HYNTY ASAKW THSGG EMFVA LNQKG IPVKG KKTKK EQKTA HFLPM AIT
MW	18.7KD
Source	E. coli cys32-thr194
Uniprot ID	P36363
Expiration Date	6 months from date of receipt.