

Product Datasheet

Guinea pig Cholecystokinin A Receptor (CCKAR) ELISA Kit orb1204423

Description	Guinea pig Cholecystokinin A Receptor (CCKAR) ELISA Kit
Reactivity	Guinea pig
Range	0.156-10ng/mL
Sensitivity	0.078ng/mL
Note	For research use only.
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q63931
Entrez	100719593
Expiration Date	Please enquire.
Application Notes	sandwich