

Product Datasheet

CD45.1 BioWave Violet 500 antibody orb1214825

Description	Mouse monoclonal antibody to CD45.1 BioWave Violet 500
Species/Host	Mouse
Reactivity	Mouse
Conjugation	violetFluor 500
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only.
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Clone ID	A20
Uniprot ID	P06800
Entrez	19264
Expiration Date	12 months from date of receipt.