

Product Datasheet

CD4 antibody (PE/Cy5) orb1214823

Description	Mouse monoclonal antibody against CD4 conjugated to PE/Cy5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only.
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Clone ID	RPA-T4
Uniprot ID	P06332
Entrez	920
Expiration Date	12 months from date of receipt.