

Product Datasheet

CD3 BioWave Violet 500 antibody orb1214821

Description	Mouse monoclonal antibody to CD3 BioWave Violet 500
Species/Host	Mouse
Reactivity	Human
Conjugation	violetFluor 500
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only.
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Clone ID	OKT3
Uniprot ID	P20963 , A8MVP6 , P09693 , P07766
Entrez	915
Expiration Date	12 months from date of receipt.