

Product Datasheet

CD8a BioWave Violet 500 antibody orb1214819

Description	Mouse monoclonal antibody to CD8a BioWave Violet 500
Species/Host	Mouse
Reactivity	Human
Conjugation	violetFluor 500
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Clone ID	OKT8
Uniprot ID	P01732
Entrez	925
Expiration Date	12 months from date of receipt.