

Product Datasheet

CD45 antibody (PerCP) orb1214818

Description	Rat monoclonal antibody against CD45 conjugated to PerCP
Species/Host	Rat
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only.
Isotype	Rat IgG2b, kappa
Clonality	Monoclonal
Clone ID	30-F11
Uniprot ID	P06800
Entrez	19264
Expiration Date	12 months from date of receipt.