

---

## Product Datasheet

### Ferret IL-8 ELISA Kit orb1216398

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Ferret IL-8 Do-It-Yourself ELISA Kit |
| <b>Reactivity</b>          | Ferret                               |
| <b>Tested Applications</b> | ELISA                                |
| <b>Target</b>              | IL-8                                 |
| <b>Storage</b>             | 4 degrees C.                         |
| <b>Note</b>                | For research use only.               |
| <b>Expiration Date</b>     | Please enquire.                      |