

Product Datasheet

SARS-CoV-2 trimeric soluble full-length Spike protein, Omicron variant, S protein orb1227191

Description	Recombinant SARS-CoV-2 trimeric soluble full-length Spike protein, Omicron variant, S protein
Reactivity	Human
Tested Applications	WB, IA, Functional Studies
Immunogen	Unconjugated
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	8xHis at the C-terminus
Note	For research use only.
Protein Sequence	16-1211
Purity	> 95%
MW	137
Source	HEK293 cells
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.