

Product Datasheet

Rabbit Anti-Duck IgG (H&L) antibody (Cy5.5) orb1289637

Description	Rabbit Polyclonal antibody against Rabbit Anti-Duck IgG conjugated to Cy5.5
Species/Host	Rabbit
Reactivity	Duck
Conjugation	Cy5.5
Tested Applications	IH, ELISA, FC, IF
Immunogen	Duck IgG
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IH: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000, E: Use at an assay dependent concentration
Expiration Date	12 months from date of receipt.