

Product Datasheet

Ferritin Light Chain antibody [Out of stock] orb1291603

Description	Ferritin Light Chain antibody
Species/Host	Rabbit
Predicted Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human Ferritin Light Chain (1-175/175 aa)
Target	Ferritin Light Chain
Form/Appearance	Liquid
Preservatives	0.01M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P02792
NCBI	2512
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.