

Product Datasheet

F(ab')₂ Fragment Goat Anti-Rabbit IgG H&L antibody (HRP) orb1291510

Description	F(ab') ₂ Fragment Goat Anti-Rabbit IgG H&L antibody (HRP)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Native Rabbit IgG
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:500-2000,IHC-P=1:400-800,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.