

Product Datasheet

WNT1 antibody orb1292558

Description	Rabbit Polyclonal antibody to WNT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	wingless-type MMTV integration site family, member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20° for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	41 kDa
Uniprot ID	P04628
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.