

Product Datasheet

NF- κ B p65 Phospho-S536 antibody orb1292536

Description	Rabbit Polyclonal antibody to NF- κ B p65 Phospho-S536
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S536 of human RelA
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20° for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	65 kDa
Uniprot ID	Q04206
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.