

Product Datasheet

ADRA2A antibody orb1292525

Description	Rabbit Polyclonal antibody to ADRA2A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Alpha-2A adrenergic receptor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20° for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	P08913
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.