

Product Datasheet

MOR antibody orb1292488

Description	Rabbit Polyclonal antibody to MOR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Mu Opioid Receptor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20° for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-70 kDa
Uniprot ID	P35372
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.