

Product Datasheet

GPR78 antibody orb1292462

Description	Rabbit Polyclonal antibody to GPR78
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	G-protein coupled receptor 78
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20° for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	39kDa
Uniprot ID	Q96P69
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.