

## Product Datasheet

### Pig C3a protein orb1291844

<b>Description</b>	Recombinant Pig C3a protein
<b>Reactivity</b>	Porcine
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	Lyophilized from a 0.2 ?m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
<b>Storage</b>	-20°C for 12 months as lyophilized;2-8°C for 1 month under sterile conditions after reconstitution
<b>Tag</b>	N-terminal His-IF2DI Tag
<b>Note</b>	For research use only.
<b>Purity</b>	Greater than 90% by SDS-PAGE gel analyses
<b>Formula</b>	Lyophilized from a 0.2 ?m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
<b>MW</b>	29.9 kDa
<b>Source</b>	E.Coli
<b>Uniprot ID</b>	<b>P01025</b>
<b>Expiration Date</b>	6 months from date of receipt.
<b>Application Notes</b>	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 ?g/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex