

Product Datasheet

Human GPRIN1 protein orb1291821

Description	Recombinant Human GPRIN1 protein
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 ?m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
Storage	-20°C for 12 months as lyophilized;2-8°C for 1 month under sterile conditions after reconstitution
Tag	N-terminal His-IF2DI Tag
Note	For research use only.
Purity	Greater than 90% by SDS-PAGE gel analyses
Formula	Lyophilized from a 0.2 ?m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
MW	45.9 kDa
Source	E.Coli
Uniprot ID	Q7Z2K8
Expiration Date	6 months from date of receipt.
Application Notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 ?g/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex