

Product Datasheet

Human LGALS2 protein orb1291816

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant Human LGALS2 protein |
| Reactivity | Human |
| Tested Applications | ELISA, WB |
| Preservatives | Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4 |
| Storage | -20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution |
| Tag | N-terminal His Tag |
| Note | For research use only. |
| Purity | Greater than 90% by SDS-PAGE gel analyses |
| Formula | Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4 |
| MW | 16.5 kDa |
| Source | E.Coli |
| Uniprot ID | P05162 |
| Expiration Date | 6 months from date of receipt. |
| Application Notes | Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex |