




Product Datasheet

FFAR2 antibody orb393286

Description	Rabbit polyclonal antibody to FFAR2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH, IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human GPR43. The exact sequence is proprietary.
Target	FFAR2
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O15552
Entrez	2867
Dilution Range	WB: 1:500:1000, IHC-P: 1:100:200, IF/ICC: 1:100:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human tonsil using FFAR2 antibody



Western blot analysis of THP1 (Lane 1) whole cell lysates using FFAR2 antibody