



Product Datasheet

CD133-1,2,3,5,7 antibody orb394857

Description	Rabbit polyclonal antibody to CD133-1,2,3,5,7.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	prominin 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20° for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	115 kDa
Uniprot ID	O43490
Dilution Range	WB: 1:1000-10000; IHC: 1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin embedded human gliomas tissue using CD133-1,2,3,5,7 antibody