

Product Datasheet

ULK2 Blocking peptide orb1440940

Description	ULK2 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only.
Protein Sequence	The synthetic peptide sequence is selected from aa 627~641 of HUMAN ULK2
Uniprot ID	Q8IYT8
NCBI	9706
Expiration Date	12 months from date of receipt.
Application Notes	Synthetic peptide