

Product Datasheet

PERK (EIF2AK3) Blocking peptide orb1440939

Description	PERK (EIF2AK3) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only.
Protein Sequence	The synthetic peptide sequence is selected from aa 482~496 of HUMAN EIF2AK3
Uniprot ID	Q9NZJ5
NCBI	9451
Expiration Date	12 months from date of receipt.
Application Notes	Synthetic peptide