

Product Datasheet

PTP1B Blocking peptide (S386) orb1440937

Description	PTP1B Blocking peptide (S386)
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only.
Protein Sequence	The synthetic peptide sequence is selected from aa 386-400 of HUMAN PTPN1
Uniprot ID	P18031
NCBI	5770
Expiration Date	12 months from date of receipt.
Application Notes	Synthetic peptide