

Product Datasheet

GSK3 β antibody orb1413523

Description	Rabbit polyclonal antibody to GSK3 β .
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	P49841 , P49840
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.