

## Product Datasheet

### Ub antibody orb1411883

<b>Description</b>	Rabbit polyclonal antibody to Ub.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Note</b>	For research use only.
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P62987</a> , <a href="#">P62979</a> , <a href="#">P0CG48</a> , <a href="#">P0CG47</a>
<b>Dilution Range</b>	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.