

Product Datasheet

NF κ B-p65 antibody orb1411432

Description	Rabbit polyclonal antibody to NF κ B-p65.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q04206
NCBI	5970
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.