

## Product Datasheet

### NCoA-3 antibody orb1411264

<b>Description</b>	Rabbit polyclonal antibody to NCoA-3.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Note</b>	For research use only.
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q9Y6Q9</a>
<b>NCBI</b>	<a href="#">8202</a>
<b>Dilution Range</b>	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.