

Product Datasheet

CD156c antibody orb1411085

Description	Rabbit polyclonal antibody to CD156c.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	O14672
NCBI	102
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.