

## Product Datasheet

### Synuclein-? antibody orb1411040

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Synuclein-?.                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                   |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                      |
| <b>Note</b>                | For research use only.                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                     |
| <b>Uniprot ID</b>          | <b>P37840</b>                                                                                                  |
| <b>NCBI</b>                | <b>6622</b>                                                                                                    |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000.<br>Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |