

Product Datasheet

IKK? antibody orb1410892

Description	Rabbit polyclonal antibody to IKK?.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q9Y6K9
NCBI	8517
Dilution Range	WB: WB 1:500-2000,IHC-p 1:500-200, ELISA 1:10000-20000
Expiration Date	12 months from date of receipt.