

Product Datasheet

Alpha-synuclein antibody orb1473902

Description	Mouse monoclonal antibody to Alpha-synuclein
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Alpha-synuclein. The exact sequence is proprietary.
Target	SNCA
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P37840 , O55042
Entrez	20617 , 6622
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.