

Product Datasheet

Human Tau protein orb1476334

| | |
|--------------------------|--|
| Description | Recombinant Human Tau protein |
| Endotoxins | Less than 1 EU/ug |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.5 |
| Preservatives | Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.5 |
| Storage | Lyophilized protein should be stored at |
| Tag | tag free |
| Note | For research use only. |
| Protein Sequence | Met1-Leu441 |
| Source | E.coli |
| Biological Origin | Human |
| Uniprot ID | P10636 |
| Expiration Date | 6 months from date of receipt. |