

Product Datasheet

ANGEL2 antibody orb1592

Description	ANGEL2 antibody
Species/Host	Rabbit
Predicted Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ANGEL2 (61-160/544 aa)
Target	ANGEL2
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	62 kDa
Uniprot ID	Q5VTE6
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:50-200, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.