

Product Datasheet

Dishevelled 3 antibody orb1327

Description	Dishevelled 3 antibody
Species/Host	Rabbit
Reactivity	Rat
Predicted Reactivity	Bovine, Canine, Gallus, Human, Mouse, Porcine, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Dishevelled 3/DVL3 (401-500/716 aa)
Target	Dishevelled 3
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	Q92997
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.