

Product Datasheet

CSFV antibody orb1388

Description	CSFV antibody
Species/Host	Rabbit
Predicted Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	whole cell protein of CSFV
Target	CSFV
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	430 kDa
Uniprot ID	P19712
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.