

Product Datasheet

IDE antibody orb10865

Description	IDE antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Gallus, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IDE (491-590/1019 aa)
Target	IDE
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	117 kDa
Uniprot ID	P14735
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.