

## Product Datasheet

### Human GPR62 Protein orb1476812

<b>Description</b>	Recombinant Human G-protein coupled receptor 62(GPR62)
<b>Endotoxins</b>	Not test.
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Preservatives</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
<b>Tag</b>	N-terminal 10xHis-tagged
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MANSTGLNASEVAGSLGLILAAVVEVGALLGNGALLVVVLRTPGLRDALY LAHLCVVDLLAAASIMPLGLLAAPPGLGRVRLGPAPCRAARFLSAALLP ACTLGVAALGLARYRLIVHPLRPGSRPPPVLVLTAVWAAAGLLGALSLLG TPPAPPPAPARCSVLGGLGPFRLWALLAFALPALLLLGAYGGIFVVAR RAALRPPRPARGSRSLHSDSLDSRLSILPPLRPRLPGGKAALAPALAVGQF AACWLPYGCACLAPAAARAAEAAAVTWVAYSFAFAHPFLYGLLQRPVRLA LGRLSRRALPGPVRACTPQAWHPRALLQCLQRPPEGPAVGPSEAPEQTPE LAGGRSPAYQGPPPESSLS-
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>MW</b>	39.1 kDa
<b>Source</b>	in vitro E.coli expression system
<b>Biological Origin</b>	Homo sapiens (Human)
<b>Uniprot ID</b>	<a href="#">Q9BZJ7</a>
<b>Expiration Date</b>	6 months from date of receipt.
<b>Application Notes</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20?/-80?. Our default final concentration of glycerol is 50%. Customers could use it as reference.