

Product Datasheet

tusA Protein orb1476737

Description	Recombinant Pseudomonas aeruginosa Sulfurtransferase TusA homolog(tusA)
Endotoxins	Not test.
Form/Appearance	Liquid or Lyophilized powder
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only.
Protein Sequence	MTHSVDAILDATGLNCPEPVMMLHNKVRDLAPGGLLKVIATDPSTRRDIP KFCVFLGHELVEQQEEAGTYLYWIRKKAD
Purity	Greater than 90% as determined by SDS-PAGE.
MW	16.3 kDa
Source	E.coli
Biological Origin	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Uniprot ID	Q9I3F3
Expiration Date	6 months from date of receipt.
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20?/-80?. Our default final concentration of glycerol is 50%. Customers could use it as reference.