

## Product Datasheet

### ISCS antibody orb819676

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| Description         | ISCS antibody                                                 |
| Species/Host        | Rabbit                                                        |
| Reactivity          | E. coli                                                       |
| Conjugation         | Unconjugated                                                  |
| Tested Applications | ELISA, WB                                                     |
| Immunogen           | Recombinant Escherichia coli O9:H4 (strain HS) iscS protein   |
| Target              | iscS                                                          |
| Form/Appearance     | Liquid                                                        |
| Preservatives       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note                | For research use only.                                        |
| Isotype             | IgG                                                           |
| Clonality           | Polyclonal                                                    |
| Uniprot ID          | <a href="#">A8A336</a>                                        |
| Expiration Date     | 12 months from date of receipt.                               |