

## Product Datasheet

### OXTR Immunizing Peptide orb1479131

|                  |                                                                  |
|------------------|------------------------------------------------------------------|
| Description      | OXTR Immunizing Peptide                                          |
| Preservatives    | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Storage          | Shipped at ambient temperature, store at -20°C                   |
| Note             | For research use only.                                           |
| Protein Sequence | ANAPKEASAFIIVC                                                   |
| Expiration Date  | 6 months from date of receipt.                                   |