

Product Datasheet

GPRC5A Antibody orb1522988

Description	GPRC5A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 307 and 357
Target	GPRC5A
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q8NFJ5
NCBI	NP_003970.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG

