

Product Datasheet

F4/80 Antibody orb1519471

Description	F4/80 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	residues 28-225 (ECD)
Target	F4/80
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Storage	2 - 8°C
Note	For research use only.
Clonality	Recombinant
Uniprot ID	Q61549
NCBI	NP_034260.1
Dilution Range	WB - 1:1,000; IP - 6 µl/1 mg lysate; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.; ICC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections.
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG