

Product Datasheet

NF-kappaB1 Antibody orb1517942

Description	NF-kappaB1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 918 and 968
Target	NF-kappaB1
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	P19838
NCBI	AAF35232.1
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG