

Product Datasheet

RelA, Phospho (S536) Antibody orb1517931

Description	RelA, Phospho (S536) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Surrounding serine 536
Target	RelA, Phospho (S536)
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q04206
NCBI	NP_068810.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG