

Product Datasheet

MADD Antibody orb1517761

Description	MADD Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1597 and 1647
Target	MADD
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q8WVG6
NCBI	NP_003673.3
Dilution Range	IHC - 1:100 - 1:500; IHC-IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG