

## Product Datasheet

### TIGAR Antibody orb1517714

|                     |                                                                      |
|---------------------|----------------------------------------------------------------------|
| Description         | TIGAR Antibody                                                       |
| Species/Host        | Rabbit                                                               |
| Reactivity          | Human                                                                |
| Conjugation         | Unconjugated                                                         |
| Tested Applications | IHC                                                                  |
| Immunogen           | Between 220 and 270                                                  |
| Target              | TIGAR                                                                |
| Form/Appearance     | Whole IgG                                                            |
| Concentration       | 250 µg/ml                                                            |
| Preservatives       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Storage             | 2 - 8°C                                                              |
| Note                | For research use only.                                               |
| Clonality           | Polyclonal                                                           |
| Uniprot ID          | <a href="#">Q9NQ88</a>                                               |
| NCBI                | <a href="#">NP_065108.1</a>                                          |
| Dilution Range      | IHC - 1:100 - 1:500                                                  |
| Expiration Date     | 12 months from date of receipt.                                      |
| Application Notes   | Format: Whole IgG                                                    |