

Product Datasheet

BAD Antibody orb1517665

Description	BAD Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1 and 50
Target	BAD
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q92934
NCBI	NP_004313.1
Dilution Range	IHC - 1:100 - 1:500; IHC-IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.
Application Notes	Format: Whole IgG