

## Product Datasheet

### CTNNBL1 Antibody orb1517554

|                     |                                                                      |
|---------------------|----------------------------------------------------------------------|
| Description         | CTNNBL1 Antibody                                                     |
| Species/Host        | Rabbit                                                               |
| Reactivity          | Human, Mouse                                                         |
| Conjugation         | Unconjugated                                                         |
| Tested Applications | IHC                                                                  |
| Immunogen           | Between 1 and 50                                                     |
| Target              | CTNNBL1                                                              |
| Form/Appearance     | Whole IgG                                                            |
| Concentration       | 250 µg/ml                                                            |
| Preservatives       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Storage             | 2 - 8°C                                                              |
| Note                | For research use only.                                               |
| Clonality           | Polyclonal                                                           |
| Uniprot ID          | <a href="#">Q8WYA6</a>                                               |
| NCBI                | <a href="#">NP_110517.2</a>                                          |
| Dilution Range      | IHC - 1:100 - 1:500                                                  |
| Expiration Date     | 12 months from date of receipt.                                      |
| Application Notes   | Format: Whole IgG                                                    |